

# GALR2 Antibody, Biotin conjugated



PACO62294

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Receptor for the hormone galanin and GALP. Receptor for the hormone spexin-1. The activity of this receptor is mediated by G proteins that activate the phospholipase C/protein kinase C pathway (via G(q)) and that inhibit adenylyl cyclase (via G(i)).

**Gene ID:**

GALR2

**Uniprot**

O43603

**Synonyms:**

Galanin receptor type 2 (GAL2-R) (GALR-2), GALR2, GALNR2

**Immunogen:**

Recombinant Human Galanin receptor type 2 protein (282-387AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A