

PPP5C Antibody, Biotin conjugated



PACO62342

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

PPP5C

Uniprot

P53041

Synonyms:

Serine/threonine-protein phosphatase 5 (PP5) (EC 3.1.3.16) (Protein phosphatase T) (PP-T) (PPT), PPP5C, PPP5

Immunogen:

Recombinant Human Serine/threonine-protein phosphatase 5 protein (370-488AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A