

PRAMEF18 Antibody, FITC conjugated



PACO62345

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

PRAMEF18

Uniprot

Q5VWM3

Synonyms:

PRAME family member 18, PRAMEF18

Immunogen:

Recombinant Human PRAME family member 18 protein (1-257AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A