

# SOX10 Antibody, HRP conjugated



PACO62380

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Transcription factor that seems to function synergistically with the POU domain protein TST-1/OCT6/SCIP. Could confer cell specificity to the function of other transcription factors in developing and mature glia.

**Gene ID:**

SOX10

**Uniprot**

P56693

**Synonyms:**

Transcription factor SOX-10, SOX10

**Immunogen:**

Peptide sequence from Human Transcription factor SOX-10 protein (49-68AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A