

AOX1 Antibody



PACO62415

Product Information

Size:

50ul

Reactivity:

Chlamydomonas reinhardtii

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000

Protein Background:

Gene ID:

AOX1

Uniprot

O65000

Synonyms:

Ubiquinol oxidase, AOX1


Immunogen:

Peptide sequence from Chlamydomonas reinhardtii Ubiquinol oxidase protein (202-217AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Western Blot. Positive WB detected in: Chlamydomonas tissue. All lanes: AOX1 antibody at 1:2000. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 39 kDa. Observed band size: 39 kDa.