

JRK Antibody



PACO62439

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:5000

Protein Background:

May bind DNA.

Gene ID:

JRK

Uniprot

O75564

Synonyms:

Jerky protein homolog, JRK, JH8


Immunogen:

Peptide sequence from Human Jerky protein homolog protein (465-483AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Western Blot. Positive WB detected in: SH-SY5Y whole cell lysate, Mouse heart tissue. All lanes: JRK antibody at 1:1000. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 62, 63 kDa. Observed band size: 62 kDa.