

RHOT2 Antibody, HRP conjugated



PACO62524

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Mitochondrial GTPase involved in mitochondrial trafficking. Probably involved in control of anterograde transport of mitochondria and their subcellular distribution.

Gene ID:

RHOT2

Uniprot

Q8IXI1

Synonyms:

Mitochondrial Rho GTPase 2 (MIRO-2) (hMiro-2) (EC 3.6.5) (Ras homolog gene family member T2), RHOT2, ARHT2 C16orf39

Immunogen:

Recombinant Human Mitochondrial Rho GTPase 2 protein (84-300AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A