

ALMS1 Antibody, HRP conjugated



PACO62564

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Involved in PCM1-dependent intracellular transport. Required, directly or indirectly, for the localization of NCAPD2 to the proximal ends of centrioles. Required for proper formation and/or maintenance of primary cilia (PC), microtubule-based structures that protrude from the surface of epithelial cells.

Gene ID:

ALMS1

Uniprot

Q8TCU4

Synonyms:

Alstrom syndrome protein 1, ALMS1, KIAA0328

Immunogen:

Recombinant Human Alstrom syndrome protein 1 protein (889-1057AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A