

# IFT74 Antibody, HRP conjugated



PACO62568

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Component of the intraflagellar transport (IFT) complex B: together with IFT81, forms a tubulin-binding module that specifically mediates transport of tubulin within the cilium. Binds beta-tubulin via its basic region. Required for ciliogenesis.

**Gene ID:**

IFT74

**Uniprot**

Q96LB3

**Synonyms:**

Intraflagellar transport protein 74 homolog (Capillary morphogenesis gene 1 protein) (CMG-1) (Coiled-coil domain-containing protein 2), IFT74, CCDC2 CMG1

**Immunogen:**

Recombinant Human Intraflagellar transport protein 74 homolog protein (375-600AA ).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A