

FCRL5 Antibody, Biotin conjugated



PACO62598

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May be involved in B-cell development and differentiation in peripheral lymphoid organs and may be useful markers of B-cell stages. May have an immunoregulatory role in marginal zone B-cells.

Gene ID:

FCRL5

Uniprot

Q96RD9

Synonyms:

Fc receptor-like protein 5 (FcR-like protein 5) (FCRL5) (BXMAS1) (Fc receptor homolog 5) (FcRH5) (Immune receptor translocation-associated protein 2) (CD antigen CD307e), FCRL5, FCRH5 IRTA2

Immunogen:

Recombinant Human Fc receptor-like protein 5 protein (419-553AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A