

VGLL1 Antibody, FITC conjugated



PACO62601

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May act as a specific coactivator for the mammalian TEFs.

Gene ID:

VGLL1

Uniprot

Q99990

Synonyms:

Transcription cofactor vestigial-like protein 1 (Vgl-1) (Protein TONDU), VGLL1, TDU

Immunogen:

Recombinant Human Transcription cofactor vestigial-like protein 1 protein (139-250AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A