

CTNNA3 Antibody, Biotin conjugated



PACO62618

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May be involved in formation of stretch-resistant cell-cell adhesion complexes.

Gene ID:

CTNNA3

Uniprot

Q9UI47

Synonyms:

Catenin alpha-3 (alpha T-catenin) (Cadherin-associated protein), CTNNA3

Immunogen:

Recombinant Human Catenin alpha-3 protein (163-302AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A