

UTP18 Antibody, HRP conjugated



PACO62660

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Involved in nucleolar processing of pre-18S ribosomal RNA.

Gene ID:

UTP18

Uniprot

Q9Y5J1

Synonyms:

U3 small nucleolar RNA-associated protein 18 homolog (WD repeat-containing protein 50), UTP18, WDR50

Immunogen:

Recombinant Human U3 small nucleolar RNA-associated protein 18 homolog protein (13-152AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A