

CAB39 Antibody, Biotin conjugated



PACO62758

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Component of a complex that binds and activates STK11/LKB1. In the complex, required to stabilize the interaction between CAB39/MO25 (CAB39/MO25alpha or CAB39L/MO25beta) and STK11/LKB1.

Gene ID:

CAB39

Uniprot

Q9Y376

Synonyms:

Calcium-binding protein 39 (MO25alpha) (Protein Mo25), CAB39, MO25

Immunogen:

Recombinant Human Calcium-binding protein 39 protein (28-290AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A