

# NSMAF Antibody, HRP conjugated



PACO62796

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Couples the p55 TNF-receptor (TNF-R55 / TNFR1) to neutral sphingomyelinase (N-SMASE). Specifically binds to the N-smase activation domain of TNF-R55. May regulate ceramide production by N-SMASE.

**Gene ID:**

NSMAF

**Uniprot**

Q92636

**Synonyms:**

Protein FAN (Factor associated with neutral sphingomyelinase activation) (Factor associated with N-SMase activation), NSMAF, FAN

**Immunogen:**

Recombinant Human Protein FAN protein (1-233AA ).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A