

# PAFAH2 Antibody, HRP conjugated



PACO62804

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Has a marked selectivity for phospholipids with short acyl chains at the sn-2 position. May share a common physiologic function with the plasma-type enzyme.

**Gene ID:**

PAFAH2

**Uniprot**

Q99487

**Synonyms:**

Platelet-activating factor acetylhydrolase 2, cytoplasmic (EC 3.1.1.47) (Serine-dependent phospholipase A2) (SD-PLA2) (hSD-PLA2), PAFAH2

**Immunogen:**

Recombinant Human Platelet-activating factor acetylhydrolase 2, cytoplasmic protein (1-206AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A