

PAFAH2 Antibody, Biotin conjugated



PACO62806

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Has a marked selectivity for phospholipids with short acyl chains at the sn-2 position. May share a common physiologic function with the plasma-type enzyme.

Gene ID:

PAFAH2

Uniprot

Q99487

Synonyms:

Platelet-activating factor acetylhydrolase 2, cytoplasmic (EC 3.1.1.47) (Serine-dependent phospholipase A2) (SD-PLA2) (hSD-PLA2), PAFAH2

Immunogen:

Recombinant Human Platelet-activating factor acetylhydrolase 2, cytoplasmic protein (1-206AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A