

ANKRD2 Antibody, FITC conjugated



PACO62817

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Functions as a negative regulator of myocyte differentiation. May interact with both sarcoplasmic structural proteins and nuclear proteins to regulate gene expression during muscle development and in response to muscle stress.

Gene ID:

ANKRD2

Uniprot

Q9GZV1

Synonyms:

Ankyrin repeat domain-containing protein 2 (Skeletal muscle ankyrin repeat protein) (hArpp), ANKRD2, ARPP

Immunogen:

Recombinant Human Ankyrin repeat domain-containing protein 2 protein (1-210AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A