

WDR35 Antibody, HRP conjugated



PACO62820

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Component of the IFT complex A (IFT-A), a complex required for retrograde ciliary transport. Required for ciliogenesis. May promote CASP3 activation and TNF-stimulated apoptosis.

Gene ID:

WDR35

Uniprot

Q9P2L0

Synonyms:

WD repeat-containing protein 35 (Intraflagellar transport protein 121 homolog), WDR35, IFT121 KIAA1336

Immunogen:

Recombinant Human WD repeat-containing protein 35 protein (954-1181AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A