

GLI2 Antibody, HRP conjugated



PACO62912

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Functions as transcription regulator in the hedgehog (Hh) pathway. Functions as transcriptional activator. May also function as transcriptional repressor. Requires STK36 for full transcriptional activator activity. Required for normal embryonic development.

Gene ID:

GLI2

Uniprot

P10070

Synonyms:

Zinc finger protein GLI2 (GLI family zinc finger protein 2) (Tax helper protein), GLI2, THP

Immunogen:

Peptide sequence from Human Zinc finger protein GLI2 protein (2-20AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A