

LRPAP1 Antibody, Biotin conjugated



PACO62942

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Molecular chaperone for LDL receptor-related proteins that may regulate their ligand binding activity along the secretory pathway.

Gene ID:

LRPAP1

Uniprot

P30533

Synonyms:

Alpha-2-macroglobulin receptor-associated protein (Alpha-2-MRAP) (Low density lipoprotein receptor-related protein-associated protein 1) (RAP), LRPAP1, A2MRAP

Immunogen:

Recombinant Human Alpha-2-macroglobulin receptor-associated protein (35-357AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A