

MYH6 Antibody, HRP conjugated



PACO62952

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Muscle contraction.

Gene ID:

MYH6

Uniprot

P13533

Synonyms:

Myosin-6 (Myosin heavy chain 6) (Myosin heavy chain, cardiac muscle alpha isoform) (MyHC-alpha), MYH6, MYHCA

Immunogen:

Recombinant Human Myosin-6 protein (894-1132AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A