

ORM1 Antibody, FITC conjugated



PACO62957

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Functions as transport protein in the blood stream. Binds various ligands in the interior of its beta-barrel domain. Also binds synthetic drugs and influences their distribution and availability in the body. Appears to function in modulating the activity of the immune system during the acute-phase reaction.

Gene ID:

ORM1

Uniprot

P02763

Synonyms:

Alpha-1-acid, glycoprotein 1 (AGP 1) (Orosomuroid-1) (OMD 1), ORM1, AGP1

Immunogen:

Recombinant Human Alpha-1-acid, glycoprotein 1 protein (19-201AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A