

SLC3A2 Antibody, Biotin conjugated



PACO62986

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

SLC3A2

Uniprot

P08195

Synonyms:

4F2 cell-surface antigen heavy chain (4F2hc) (4F2 heavy chain antigen) (Lymphocyte activation antigen 4F2 large subunit) (Solute carrier family 3 member 2) (CD antigen CD98), SLC3A2, MDU1

Immunogen:

Recombinant Human 4F2 cell-surface antigen heavy chain protein (404-630AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A