

# SRPX Antibody, Biotin conjugated



PACO62994

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

May be involved in phagocytosis during disk shedding, cell adhesion to cells other than the pigment epithelium or signal transduction.

**Gene ID:**

SRPX

**Uniprot**

P78539

**Synonyms:**

Sushi repeat-containing protein SRPX, SRPX, ETX1

**Immunogen:**

Recombinant Human Sushi repeat-containing protein SRPX protein (31-237AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A