

# benC Antibody



PACO63027

---

## Product Information

**Size:**

50ul

**Reactivity:**

Streptomyces sp

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:1000-1:5000

**Protein Background:**

Carrier of the growing fatty acid, chain in fatty acid, biosynthesis.

**Gene ID:**

benC

**Uniprot**

B1GSN3

**Synonyms:**

BenC protein, benC

**Immunogen:**


Recombinant Streptomyces sp BenC protein (1-86AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---



Western Blot. Positive WB detected in Recombinant protein. All lanes: benC antibody at 1:2000. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 33 kDa. Observed band size: 33 kDa.