

NKAPL Antibody, FITC conjugated



PACO63101

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Transcriptional repressor of Notch-mediated signaling. Required for spermatogenesis.

Gene ID:

NKAPL

Uniprot

Q5M9Q1

Synonyms:

NKAP-like protein, NKAPL, C6orf194

Immunogen:

Peptide sequence from Human NKAP-like protein (114-123AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A