

# CTHRC1 Antibody, FITC conjugated



PACO63161

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

May act as a negative regulator of collagen matrix deposition.

**Gene ID:**

CTHRC1

**Uniprot**

Q96CG8

**Synonyms:**

Collagen triple helix repeat-containing protein 1 (Protein NMTC1), CTHRC1

**Immunogen:**

Recombinant Human Collagen triple helix repeat-containing protein 1 protein (31-243AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A