

# TMEM237 Antibody, HRP conjugated



PACO63168

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Component of the transition zone in primary cilia. Required for ciliogenesis.

**Gene ID:**

TMEM237

**Uniprot**

Q96Q45

**Synonyms:**

Transmembrane protein 237 (Amyotrophic lateral sclerosis 2 chromosomal region candidate gene 4 protein), TMEM237, ALS2CR4

**Immunogen:**

Peptide sequence from Human Transmembrane protein 237 protein (59-78AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A