

MIOS Antibody, Biotin conjugated



PACO63210

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

As a component of the GATOR subcomplex GATOR2, functions within the amino acid, sensing branch of the TORC1 signaling pathway. Indirectly activates mTORC1 and the TORC1 signaling pathway through the inhibition of the GATOR1 subcomplex. It is negatively regulated by the upstream amino acid, sensors SESN2 and CASTOR1.

Gene ID:

MIOS

Uniprot

Q9NXC5

Synonyms:

GATOR complex protein MIOS (Missing oocyte meiosis regulator homolog), MIOS

Immunogen:

Recombinant Human GATOR complex protein MIOS protein (1-237AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A