

TJP2 Antibody, HRP conjugated



PACO63232

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Plays a role in tight junctions and adherens junctions.

Gene ID:

TJP2

Uniprot

Q9UDY2

Synonyms:

Tight junction protein ZO-2 (Tight junction protein 2) (Zona occludens protein 2) (Zonula occludens protein 2), TJP2, X104 ZO2

Immunogen:

Peptide sequence from Human Tight junction protein ZO-2 protein (1165-1183AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A