

SH2D3A Antibody, HRP conjugated



PACO63236

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May play a role in JNK activation.

Gene ID:

SH2D3A

Uniprot

Q9BRG2

Synonyms:

SH2 domain-containing protein 3A (Novel SH2-containing protein 1), SH2D3A, NSP1

Immunogen:

Peptide sequence from Human SH2 domain-containing protein 3A protein (211-231AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A