

WWC3 Antibody, FITC conjugated



PACO63249

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

cytosol, kinase binding, protein complex scaffold, negative regulation of hippo signaling, negative regulation of organ growth, negative regulation of transcription from RNA polymerase II promoter.

Gene ID:

WWC3

Uniprot

Q9ULE0

Synonyms:

Protein WWC3, WWC3, KIAA1280

Immunogen:

Peptide sequence from Human Protein WWC3 protein (65-85AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A