

CRTAC1 Antibody, Biotin conjugated



PACO63282

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

extracellular exosome.

Gene ID:

CRTAC1

Uniprot

Q9NQ79

Synonyms:

Cartilage acid, c protein 1 (68 kDa chondrocyte-expressed protein) (CEP-68) (ASPIC), CRTAC1, ASPIC1 CEP68

Immunogen:

Recombinant Human Cartilage acid, c protein 1 protein (309-661AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A