

# ERC1 Antibody, FITC conjugated



PACO63301

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Regulatory subunit of the IKK complex. Probably recruits I $\kappa$ B $\alpha$ /NF $\kappa$ BIA to the complex. May be involved in the organization of the cytomatrix at the nerve terminals active zone (CAZ) which regulates neurotransmitter release. May be involved in vesicle trafficking at the CAZ. May be involved in Rab-6 regulated endosomes to Golgi transport.

**Gene ID:**

ERC1

**Uniprot**

Q8IUD2

**Synonyms:**

ELKS/Rab6-interacting/CAST family member 1 (ERC-1) (Rab6-interacting protein 2), ERC1, ELKS KIAA1081 RAB6IP2

**Immunogen:**

Recombinant Human ELKS/Rab6-interacting/CAST family member 1 protein (1-224AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A