

ERC1 Antibody, Biotin conjugated



PACO63302

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Regulatory subunit of the IKK complex. Probably recruits I κ B α /NF κ BIA to the complex. May be involved in the organization of the cytomatrix at the nerve terminals active zone (CAZ) which regulates neurotransmitter release. May be involved in vesicle trafficking at the CAZ. May be involved in Rab-6 regulated endosomes to Golgi transport.

Gene ID:

ERC1

Uniprot

Q8IUD2

Synonyms:

ELKS/Rab6-interacting/CAST family member 1 (ERC-1) (Rab6-interacting protein 2), ERC1, ELKS KIAA1081 RAB6IP2

Immunogen:

Recombinant Human ELKS/Rab6-interacting/CAST family member 1 protein (1-224AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A