

LPAR2 Antibody, HRP conjugated



PACO63344

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Receptor for lysophosphatidic acid, (LPA), a mediator of diverse cellular activities. Seems to be coupled to the G(i)/G(o), G(12)/G(13), and G(q) families of heteromeric G proteins. Plays a key role in phospholipase C-beta (PLC-beta) signaling pathway. Stimulates phospholipase C (PLC) activity in a manner that is independent of RALA activation.

Gene ID:

LPAR2

Uniprot

Q9HBW0

Synonyms:

Lysophosphatidic acid, receptor 2 (LPA receptor 2) (LPA-2) (Lysophosphatidic acid, receptor Edg-4), LPAR2, EDG4 LPA2

Immunogen:

Peptide sequence from Human Lysophosphatidic acid, receptor 2 protein (8-27AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A