

RHBDD3 Antibody, HRP conjugated



PACO63404

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

RHBDD3

Uniprot

Q9Y3P4

Synonyms:

Rhomboid domain-containing protein 3, RHBDD3, C22orf3

Immunogen:

Peptide sequence from Human Rhomboid domain-containing protein 3 protein (304-323AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A