

# APCDD1 Antibody, HRP conjugated



PACO63444

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Negative regulator of the Wnt signaling pathway. Inhibits Wnt signaling in a cell-autonomous manner and functions upstream of beta-catenin. May act via its interaction with Wnt and LRP proteins. May play a role in colorectal tumorigenesis.

**Gene ID:**

APCDD1

**Uniprot**

Q8J025

**Synonyms:**

Protein APCDD1 (Adenomatosis polyposis coli down-regulated 1 protein), APCDD1, DRAPC1

**Immunogen:**

Recombinant Human Protein APCDD1 protein (22-288AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A