

# OSBPL7 Antibody, HRP conjugated



PACO63488

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

autophagosome, cytosol, perinuclear endoplasmic reticulum, plasma membrane, cholesterol binding, sterol transporter activity, bile acid, biosynthetic process, cellular response to cholesterol, positive regulation of proteasomal protein catabolic process, regulation of autophagy.

**Gene ID:**

OSBPL7

**Uniprot**

Q9BZF2

**Synonyms:**

Oxysterol-binding protein-related protein 7 (ORP-7) (OSBP-related protein 7), OSBPL7, ORP7

**Immunogen:**

Recombinant Human Oxysterol-binding protein-related protein 7 protein (1-249AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A