

SNRK Antibody, HRP conjugated



PACO63492

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May play a role in hematopoietic cell proliferation or differentiation. Potential mediator of neuronal apoptosis.

Gene ID:

SNRK

Uniprot

Q9NRH2

Synonyms:

SNF-related serine/threonine-protein kinase (EC 2.7.11.1) (SNF1-related kinase), SNRK, KIAA0096 SNFRK

Immunogen:

Peptide sequence from Human SNF-related serine/threonine-protein kinase protein (340-360AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A