

TENM1 Antibody, Biotin conjugated



PACO63526

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Involved in neural development, regulating the establishment of proper connectivity within the nervous system. May function as a cellular signal transducer.

Gene ID:

TENM1

Uniprot

Q9UKZ4

Synonyms:

Teneurin-1 (Ten-1) (Protein Odd Oz/ten-m homolog 1) (Tenascin-M1) (Ten-m1) (Teneurin transmembrane protein 1) [Cleaved into: Ten-1 intracellular domain (IDten-1) (Ten-1 ICD); Teneurin C-terminal-associated peptide (TCPA-1) (Ten-1 extracellular domain) (Ten-1 ECD)], TENM1, ODZ1 TNM1

Immunogen:

Recombinant Human Teneurin-1 protein (1-241AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A