

ARL2BP Antibody, HRP conjugated



PACO63536

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Together with ARL2, plays a role in the nuclear translocation, retention and transcriptional activity of STAT3. May play a role as an effector of ARL2.

Gene ID:

ARL2BP

Uniprot

Q9Y2Y0

Synonyms:

ADP-ribosylation factor-like protein 2-binding protein (ARF-like 2-binding protein) (ARL2-binding protein) (Binder of ARF2 protein 1), ARL2BP, BART BART1

Immunogen:

Recombinant Human ADP-ribosylation factor-like protein 2-binding protein (1-163AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A