

# ifnlr1 Antibody, HRP conjugated



PACO63560

---

## Product Information

**Size:**

50ul

**Reactivity:**

Zebrafish

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

cytokine-mediated signaling pathway.

**Gene ID:**

ifnlr1

**Uniprot**

D9IL84

**Synonyms:**

Cytokine receptor family member b14, ifnlr1, crfb14

**Immunogen:**

Peptide sequence from Zebrafish Cytokine receptor family member b14 protein (106-130AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A