

ifnlr1 Antibody, Biotin conjugated



PACO63562

Product Information

Size:

50ul

Reactivity:

Zebrafish

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

cytokine-mediated signaling pathway.

Gene ID:

ifnlr1

Uniprot

D9IL84

Synonyms:

Cytokine receptor family member b14, ifnlr1, crfb14

Immunogen:

Peptide sequence from Zebrafish Cytokine receptor family member b14 protein (106-130AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A