

# CD4L-2 Antibody, HRP conjugated



PACO63568

---

## Product Information

**Size:**

50ul

**Reactivity:**

Grass carp

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

CD4L-2

**Uniprot**

A0A1B2FIB0

**Synonyms:**

CD4L-2

**Immunogen:**

Peptide sequence from Grass carp CD4-like protein 2 protein (59-77AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A