

CD4L-2 Antibody, Biotin conjugated



PACO63570

Product Information

Size:

50ul

Reactivity:

Grass carp

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

CD4L-2

Uniprot

A0A1B2FIB0

Synonyms:

CD4L-2

Immunogen:

Peptide sequence from Grass carp CD4-like protein 2 protein (59-77AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A