

plac8.1 Antibody



PACO63571

Product Information

Size:

50ul

Reactivity:

Zebrafish

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

plac8.1

Uniprot

Q498P1

Synonyms:

Placenta-specific 8, tandem duplicate 1, Zgc:114201, plac8.1, zgc:114201

Immunogen:

Peptide sequence from Zebrafish Placenta-specific 8, tandem duplicate 1 protein (66-82AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A