

plac8.2 Antibody



PACO63575

Product Information

Size:

50ul

Reactivity:

Zebrafish

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

plac8.2

Uniprot

A2VD52

Synonyms:

Zgc:158845 protein, plac8.2

Immunogen:

Peptide sequence from Zebrafish Zgc:158845 protein (5-19AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A