

CAMK2N1 Antibody, HRP conjugated



PACO63628

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Potent and specific inhibitor of CaM-kinase II (CAMK2).

Gene ID:

CAMK2N1

Uniprot

Q7Z7J9

Synonyms:

Calcium/calmodulin-dependent protein kinase II inhibitor 1 (CaMKII inhibitory protein alpha) (CaMKIIN-alpha), CAMK2N1

Immunogen:

Peptide sequence from Human Calcium/calmodulin-dependent protein kinase II inhibitor 1 protein (9-27AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A